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This surgical technique is exclusively intended for medical professionals, especially physicians and surgeons.
This document does not constitute medical advice, it does not dispense medical recommendations and it does not convey any diagnostic or therapeutic information.
Informations and techniques presented in this document were compiled by a team of medical experts and Gruppo Bioimpianti's specialists;
however Gruppo Bioimpianti excludes any liability for improper use of informations.

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Müller type cups feature hemispherical shape manufactured entirely from ultra-high-molecular-weight polyethylene (UHMWPE).

They are cemented cups, which means that on their convex external surface they have a special grooved design aimed at improving the cup's primary fastening. This primary fastening takes place on the cup's impact on the acetabular cavity, which is previously filled with bone cement. The purpose of these grooves is also to increase the specific surface on contact with the acetabulum, providing greater peripheral contact and limiting the appearance of possible micro-movements that may cause a mobilization.

The hemispherical design enables the maximum bone stock to be conserved, allowing uniform load distribution to the bone thanks to a narrow contact with the acetabular bone.

All the cups have a marker ring whose purpose is to assist the surgeon in checking the correct acetabular placement of the implanted cup. Since UHMWPE is a radiotransparent material, this opaque-radio marker ring is necessary for later post-operative clinical monitoring.

The external holes (three on the edge of Müller Standard, Snap-fit and Double Snap-fit cups, one on Müller with 10° Antiluxation shoulder) are for lodging the specific instruments to facilitate positioning, orientation and impaction on the acetabulum.

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Adattatore 11-06100xx
Adaptor 11-06100xx
Adaptador 11-06100xx

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Place Hohmann separators on the anterior and posterior edges and below the acetabular groove in order to obtain a complete view of the cavity.

Use the **Acetabolar Broach** to prepare the acetabulum by eliminating the acetabular cartilage until there is a bleeding subchondral bone. (Fig. 1).

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After meticulously washing and cleaning the acetabulum, insert a portion of the acrylic cement in a paste state at pressure into the cavity.

The Müller cup is mounted in **Müller Cup Impactor 11-0610000** and the **Adaptor 11-06100xx** (Fig. 2-3) that have to be chosen depending on the type of cup (Charts page 5).

Press against the bone wall of the acetabulum to express the excess cement. Keep the **Müller Cup Impactor** with the acetabular cup immobile until the cement has set. The implant selected must have a diameter of 2 or 4 mm less than the acetabular cavity in order to obtain a sufficient thickness of cement layer.

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REDUCTION

Remove the exceeding cement from the edges of the acetabulum and eliminate fragments from the interior of the cup, before reducing the femoral head.

When the cement is completely hard, proceed to perform a reduction with the implanted stem and a trial head.

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ONLY FOR THE SNAP-FIT CUP

To reduce the Müller Snap Fit Cup, the cup and the stem must be implanted before inserting the femoral head on the stem.

Press the femoral head (placed on the neck of the stem) into the cup until it snaps in place. Make sure that the femoral head is fully inserted into the acetabulum and that there is no difficulty in movement.

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ONLY FOR THE DOUBLE SNAP-FIT CUP

To reduce the Müller Double Snap Fit Cup, the cup and the stem must be implanted before inserting the femoral head on the stem.

Press the femoral head (placed on the neck of the stem) into the cup until it snaps in place. Make sure that the femoral head is fully inserted into the acetabulum and that there is no difficulty in movement.

Finally insert the reinforcement ring as shown in figure 4.

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REF: 11-0610000

2 ADATTATORE/ ADAPTOR/ ADAPTADOR

REF.	DESCRIPTION
11-0610010	Adaptor for Muller Cup 28mm
11-0610020	Adaptor for Muller Cup 32mm
11-0610025	Adaptor for Muller Cup Low Profile 28mm
11-0610030	Adaptor for Muller Cup 28mm 10° 42-50
11-0610035	Adaptor for Muller Cup 28mm 10° 52-60
11-0610040	Adaptor for Muller Cup 32mm 10° 44-50
11-0610045	Adaptor for Muller Cup 32mm 10° 52-60
11-0610050	Adaptor for Muller Snap Fit 28mm

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All components of the implant must be handled and stored with extreme care, as any alteration (deformation, scratches, and so on) could put the implant at risk. In the operating theatre it is likewise necessary to make sure the cup isn't impaired by the instruments. For reasons of compatibility between materials and geometrical form, only GRUPPO BIOIMPIANTI instruments should be used. As with implant materials, great care must also be taken with the instruments, which must be stored in their racks to prevent accidental impact and surface damage.

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Each cup is supplied appropriately packaged and marked ensuring error-free identification and sterilized with Ethylene Oxide (EtO).

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