# SMR 100's

MODULAR HIP STEM  
VASTAGO MODULAR DE CADERA

OPERATING TECHNIQUE

TÉCNICA QUIRÚRGICA

EN | ES

*Enjoy Mobility*

# SMR 100's

MODULAR HIP STEM  
VASTAGO MODULAR DE CADERA

## SMR 100'S

The SMR 100'S stem is a modular system designed for hip replacement in cases of dysplasia, results of fractures and osteotomies. The stem features three components: proximal component, distal component and locking screw. The proximal component is available in sizes from 40mm to 70mm, The distal component is 100mm long.

## Indications

The SMR 100's stem is indicated for use in partial or total hip arthroplasty for the treatment of the following diseases: non-inflammatory degenerative joint disease such as osteoarthritis, avascular necrosis and hip dysplasia; rheumatoid arthritis; osteoarthritis after femoral heads fractures; results of osteotomy; femoral canals with particular morphologies. Furthermore, this stem can be used for the correction of functional deformities.

This stem is indicated for cementless implants.

## Contraindications

The hip joint surgery is absolutely contraindicated in cases of: systemic or local infection, sepsis, and osteomyelitis.

It is relatively contraindicated in case of: osteoporosis; patient uncooperative or suffering from neurological disorders, unable to follow directions; systemic disorders and / or metabolic problems that lead to a progressive deterioration of bone support; neurological or neuromuscular disorders that could create an unacceptable risk to the prostheses instability or lead to a failure of prostheses fixation; osteomalacia; active infection or suspected latent infection in the hip joint; distant focus of infection that could spread to the implant site; persistent acute or chronic osteomyelitis; vascular insufficiency, muscular atrophy, neuromuscular diseases; incomplete or insufficient presence of soft tissue around the knee joint; obesity; inadequate bone stock for the prostheses support or fixation; skeletal immaturity local or disseminated neoplastic diseases; incorrigible severe deformities.

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## PROXIMAL COMPONENT

EN

The modular proximal component solves the problem of restoring the physiological length of the limb and stabilizing the joint thanks to its four lengths 40, 50, 60 and 70mm and two CCD angles 125° and 135°.

The mirror polished neck area has been studied to decrease polyethylene wear during flexion. In addition, the mirror finish eliminates the risk of fatigue cracks.

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## DISTAL COMPONENT

EN

The distal component is conical and has longitudinal fins, this prosthesis design guarantees axial stability and resistance to torsional forces at the same time.

The modular design has nine diameters (14, 15, 16, 17, 18, 19, 20, 22 and 24mm).

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## ASSEMBLY

EN

The two components are assembled using a "Morse taper" 2°51' so that the proximal component can be turned into position with the anteversion required by the surgeon. The system delivers continuous 360° rotation.

The locking screw is used to hold the "Morse taper" cone in place but does not have a sealing effect, as it is the coupling of the two tapers that guarantees the seal over time. The screw has a cross bolt made of polyethylene to prevent loosening.

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## MATERIALS

EN

Titanium alloy Ti 6Al 4V Grade 5 ELI, which ensures excellent biocompatibility and mechanical strength. The surface finish is obtained by abrasion with corundum to increase the contact surface and thereby improve prosthesis-bone anchorage.

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### Pre-operative assessment

EN

Thorough pre-operative planning enables the surgeon to adopt a correct approach to the operation and foresee certain situations that might arise.

The objective is to choose the diameter and length of the distal component and size of the proximal component. Pre-operative choice of components is indicative, as the definitive size will be determined intra-operatively.

Pre-operative planning requires an anterior-posterior X-ray of the pelvis and pre-operative templates of the SMR 100's stem (provided in the instrument kit). The templates are 15% enlargements. (FIG.1)

Align the template with the anatomical axis of the femur to determine the size of the distal component, making sure that it's well anchored to good quality cortical bone.

Any anatomical differences may be corrected using the femoral heads and proximal components, which come in different lengths. Also determine the exact position of the longitudinal osteotomy, where possible.

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### Reaming of distal component

EN

Having completed the removal of the previous implant, proceed to reaming with the conical reamers (Ref. 120540714 ÷ 120540724\*). (FIG.2)

Reaming may be done manually with the "T" handle (Ref. 110381521) or with a motor, using the Jacobs Adaptor (Ref. 110381525).

The instrument kit has reamers with diameters from 14 mm to 20 mm, in 1 mm steps, and 22 mm and 24 mm.

Start with a reamer that's one or two millimetres smaller than the measurement taken in pre-operative planning.

The reamers must be inserted into the femoral canal till the depth indicator is aligned with the apex of the greater trochanter.

The shaft of the reamers is graduated and divided from 40 to 70, corresponding to the four available sizes of the proximal components.

For example, if the reamer stops at 60, a 60mm proximal component will be used. Work with conical reamers of successive sizes till achieving optimum stability and resistance and coming into contact with the cortical.

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## 3 Insertion of distal component

EN

Use of the reamers yields information on the diameter and correct depth of the distal component in the cortical.

Mount the distal component inserter (Ref. 120540410) on the chosen component (Ref. 120373914 ÷ 120373924\*) and proceed to feed it into the previously prepared canal.

Proceed to insertion with the distal component inserter advancing gradually and checking the depth, as far as the maximum, to obtain good torsional stability and avoiding any unexpected penetration under load.

(FIG.3)

Indicatively, the final position of the component will correspond to the mark identified on the reamer in the previous phase.

Any difference may be corrected using the proximal components and femoral heads, which come in different lengths.

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
FIG.3

## 4 Preparation of the bone for the proximal component

EN

Remove the distal component inserter (Ref. 120540410) from the distal component, screw the guide shaft for tubular reamer (Ref. 120540415) onto the distal component. (FIG.4a)

When carrying out this operation, make sure it's done in the right direction, as a mistake could irreparably damage the distal component and compromise all the successive phases.

Now ream the proximal portion with the appropriate tubular reamer (Ref. 120540400) to prepare a housing for the proximal component. (FIG.4b)

Attention: proceed with prudence in the bone reaming phase to avoid damaging any soft tissue in the field of action.

Attention: the reamer must rest on the apex of the Guide in order to avoid difficulties with insertion in the subsequent housing of the proximal components. To be certain of carrying out the procedure correctly, make sure there is no residue inside the reamer that might hinder its correct penetration in the bone. There is an inspection hole on the reamers for checking there are no residues.


FIG.4a

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# SMR 100's

MODULAR HIP STEM  
VASTAGO MODULAR DE CADERA

## 5 Insertion of trial proximal component

EN

Insert the trial proximal component and lock it with the appropriate screw (Ref. 120540040 ÷ 120540070). Do not tighten excessively. (FIG.5a)

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Carry out the trial reduction and assess the stability of the hip, the length of the limb and the breadth of movement. Also check for the desired anteversion, marking the bone by the reference line marked laterally on the trial proximal component.

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FIG. 5a

FIG. 5b

## 7 Insertion of the definitive proximal component

EN

Remove the trial head and trial proximal component using the revision proximal component extractor (Ref. 120540420).

Proceed with prudence in the extraction phase and make sure the correct direction is chosen, as a mistake could compromise all the successive phases.

Take the definitive proximal component (having the same size as the trial one; Ref. 120362540 ÷ 120362570; 120363540 ÷ 120363570) and use the revision proximal component impactor (Ref. 120540411) to insert it in the Morse taper of the previously cleaned distal component.

(FIG.6)

Be careful to observe the correct anteversion by aligning the reference line on the proximal component with the marking previously made on the bone.

Definitively impact the proximal component.

## Insertion of the definitive proximal component

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## 8 Blocking the proximal component

EN

Insert the locking screw (Ref. 120376140 ÷ 120376170) of the appropriate size for the proximal component used and tighten it with the screwdriver for revision/resection prosthesis (Ref. 120540405).

When tightening the screw, use the non-rotation device (Ref. 120540412), which fits into the Morse taper, to prevent rotation during assembly.

(FIG.7)

The locking screw ensures additional anchoring between the distal and proximal components, as well as holding the Morse taper.

Insert the definitive head (corresponding to the trial used) by impacting it with the appropriate adaptor for head impactor 36mm (Ref. 110380812) mounted on the head impactor (Ref. 110380800) and reduce the prosthesis.

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
FIG.7

## 9 Extraction of SMR 100's proximal component

EN

Should it be necessary to extract the proximal component of the SMR 100's stem, remove the locking screw, screw the extractor for revision proximal component (Ref. 120540420) into the locking screw hole and turn until it's fully extracted. (FIG.8)

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## 10 Extraction of SMR 100's distal component

EN

Should it be necessary to extract the distal component of the SMR 100's Stem, use the instruments provided for the purpose. Screw the extension for distal component extraction ( Ref. 120540421) onto the distal component and use the slap hammer (Ref. 140030665) to extract it.

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## Warnings EN

Never combine components made by different manufacturers because there may be different tolerances that might prejudice the outcome of the implant. In particular, never use stems and femoral heads produced by different manufacturers, because coupling tapers can be non-compatible. Gruppo Bioimpianti S.r.l. will not respond for any problems arising by a similar choice. The surgeon must be aware that a damaged or scratched neck can influence stem resistance and lead to a premature neck fracture.

1. Handle the product with due care to avoid accidental scratches, even during surgery. An inadequate handling of the device, pre-operative and / or intra-operative, or any damage can lead to corrosion, fatigue failure or excessive wear. Use clean gloves when handling the device.
2. Use screws or accessories made of suitable material - in the case of a system composed of several elements, always use the same kit elements, in case of acetabular cups made of titanium alloy (ISO 5832/3) use, if necessary, only screws of the same material.
3. Follow the surgical technique provided by Gruppo Bioimpianti.
4. A complete cleaning and removal of metal particles, bone cement, as well as any other debris resulting from surgery, it is necessary to minimize the wear of the articular surface.
5. Do not use the device if there is any damage.
6. Do not alter the system in any way.
7. Inadequate choice, positioning, alignment and fixation of one or more components could result in an incorrect distribution of the stresses, with a consequent reduction of life of the prosthetic implant.
8. An incorrect alignment or a plant not accurate can lead to a phenomenon of excessive wear and / or failure of the system itself.
9. All implantable devices for hip replacement are DISPOSABLE so it's not allowed their reuse. The repeated use leads to a system burdened by a previous use, which affects the fatigue strength and thus increases the risk of breakage.
10. Not implant a component in a patient if it had already been implanted, albeit temporarily, in another patient.
11. Utilize siempre el tornillo de seguridad adecuado en el vástago SMR (revisión) y MRP (tumor).
12. Para el vástago SMR, el tornillo de seguridad debe tener la misma longitud del componente proximal.
13. It is the responsibility of the operating surgeon to determine whether there is adequate initial fixation and stability
14. Patient smoking may result in a delayed healing, non-healing and/or compromised stability.
15. It is the responsibility of the operating surgeon to ensure that patient is not sensitive to the materials used to manufacture the implant's components.

Implantable devices for hip replacement surgery provides the surgeons the necessary tools to reduce pain and restore the Hip joint function. While this product is generally successful in attaining these goals, it cannot be expected to withstand the activity levels and function of a normal healthy joint.

## Precautions EN

The instruments set, realized by Gruppo Bioimpianti, was designed and manufactured to provide surgeons with the necessary instruments to realize a simple and safe implantation. The use of components and/or instruments not included in the instruments set provided by Gruppo Bioimpianti S.r.l., except the generic instruments present in the operating rooms (drill, oscillating blades, ) could determine negative events. Any fracture or breaking of instruments shall be reported. Surgical instruments are inevitably subject to wear with normal usage and if they have experienced excessive force are susceptible to fracture. Surgical instruments should only be used for their intended purpose, as illustrated in the surgical technique provided by Gruppo Bioimpianti. It is recommended that all instruments be regularly inspected to identify any possible damage. Do not leave in the surgical site any trial, packaging, instrument or its part.

## Possible adverse effects EN

The possible adverse effects that may result from the use of a prosthetic system, are the classics that may result from a hip arthroplasty and include the following:

1. Osteolysis (progressive bone resorption). Osteolysis can be asymptomatic and therefore routine periodic radiographic examination is essential to prevent any serious future complication
2. Particles leading to increased wear rates necessitating early revision.
3. Early or late deep infection which may necessitate the removal of the prosthesis.
4. Allergic reactions.
5. Intra-operative bone fractures particularly in the presence of low bone density due to osteoporosis, bone defects resulting from previous surgeries, bone resorption, or insertion of the device.
6. Damage to blood vessels or hematoma
7. Dislocation, migration and / or subluxation of the prosthetic components due to trauma, incorrect attachment, loss of fixation, misalignment, improper positioning, bone resorption, soft tissue lacyden, sudden movement and / or natural and / or excessive activity.
8. Peri-articular calcification or ossification, with or without impediment to joint mobility.
9. Undesirable variation of the length of the limb.
10. Reduced Range of Movement due to an incorrect choice of the components or a inadequate placement of the implant, or femoral impingement.
11. Fatigue fracture of the prosthetic components can occur as result of loss of fixation or incomplete fixation, strenuous activity, trauma, misalignment or excessive weight.
12. Corrosion at the interface between components.
13. Wear and / or deformation of the articular surfaces.
14. Pain.
15. Post-operative bone fracture post-operative and / or post-operative pain.
16. Audible sounds during motion.
17. Nerve damage.
18. Trochanteric avulsion or non-union as result of excess muscular tension, early weight bearing, or inadequate reattachment.
19. Problems of the knee or ankle of the affected limb or contralateral limb aggravated by leg length discrepancy, too much femoral medialization or muscle deficiencies.

## Sterilization EN

Product is supplied sterile (according to the definition reported on the standard UNI EN 556). In particular metal components are sterilized through rays (Beta or Gamma) according to the normative UNI EN ISO 11137-1,-2,-3 while the elements manufactured in UHMWPE are sterilized with ethylene oxide (ETO) according the normative UNI EN ISO 11135-1. Both sterilizations are guaranteed for 10 years. The sterilization method is indicated on the label in which is also reported the expiry date that indicates the maximum time to use the product. The packaging (double pack) allows products to be handled by non-sterile personnel during conservation. The extraction of the product must be carried out by sterile operating personnel. If at any time during unpacking an anomalous situation arises, it is advisable not to use the piece because there is a risk that it may have lost the sterility and its return to Gruppo Bioimpianti which will carry out, if it is possible, a reconditioner.

## Identification and traceability EN

All products manufactured by Gruppo Bioimpianti s.r.l. are identified by the label on the packaging which states the contents, size or another parameter which indicates the dimensions, the product's code, the production batch, the sterility symbol with the expiration date and the sterilization method. To guarantee external traceability even after implantation they are provided six labels inside the packaging, which shall be used on the surgery register and on the patients' medical records.

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## Disclaimer EN

This surgical technique is exclusively intended for medical professionals, especially physicians and surgeons. This document does not constitute medical advice, it does not dispense medical recommendations and it does not convey any diagnostic or therapeutic information.

Informations and techniques presented in this document were compiled by a team of medical experts and Gruppo Bioimpianti's specialists, however Gruppo Bioimpianti excludes any liability for improper use of informations.

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