



Quadrofit Total Knee System




Your partner for implants and
trauma products

Quadrofit Total Knee System

The Quadrofit Total Knee System is engineered for both cemented and cementless total knee arthroplasty. Utilizing only one surgical instrument set, it allows surgeons to intraoperatively mix and match components to meet a wide range of patient needs.

The 4 Surgical Combinations

	<i>Femoral Component</i>	<i>Tibial Component</i>	<i>Tibial Insert</i>
<i>Option 1:</i>	<i>Standard</i>	<i>Standard</i>	<i>Standard Fixed</i>
<i>Option 2:</i>	<i>Standard</i>	<i>Mobile Bearing (MB)</i>	<i>Standard MB</i>
<i>Option 3:</i>	<i>PS</i>	<i>Standard</i>	<i>PS Fixed</i>
<i>Option 4:</i>	<i>PS</i>	<i>Mobile Bearing (MB)</i>	<i>PS MB</i>



Option 1: Standard Fixed Option 2: Standard MB Option 3: PS Fixed Option 4: PS MB

Key Clinical Features

- **Bone-Preserving Femoral Design:** Minimal resections adapted to different options to reduce bone loss
- **Anatomical Design:** An optimized trochlea reduces tendon pressure and helps prevent patellar dislocation
- **Revision Solutions:** Tibial Component Half Spacers and Tibia stems are available to compensate for bone loss

Material Specifications

- **Femoral & Tibial Components:** CoCrMo alloy (ASTM F75 - ISO 5832/4)
- **Tibial Inserts & Patella:** PE-UHMW-Cross Linked and PE-UHMW (ISO 5834/1-2)

Femoral Component



Quadrofit Standard Femoral Component (Sterile)				
Cemented CoCr			Cementless CoCr Ti + Ha Coated	
Right	Left	Size	Right	Left
EB8101C000S	EB8102C000S	0	EB8131C000S	EB8132C000S
EB8101C001S	EB8102C001S	1	EB8131C001S	EB8132C001S
EB8101C002S	EB8102C002S	2	EB8131C002S	EB8132C002S
EB8101C003S	EB8102C003S	3	EB8131C003S	EB8132C003S
EB8101C004S	EB8102C004S	4	EB8131C004S	EB8132C004S
EB8101C005S	EB8102C005S	5	EB8131C005S	EB8132C005S

Quadrofit Posterior Stabilized Femoral Component (Sterile)					
Cemented CoCr			Cementless CoCr Ti + Ha Coated		
Right	Left	Size	Right	Left	
EB8121C000S	EB8122C000S	0	EB8135C000S	EB8136C000S	
EB8121C001S	EB8122C001S	1	EB8135C001S	EB8136C001S	
EB8121C002S	EB8122C002S	2	EB8135C002S	EB8136C002S	
EB8121C003S	EB8122C003S	3	EB8135C003S	EB8136C003S	
EB8121C004S	EB8122C004S	4	EB8135C004S	EB8136C004S	
EB8121C005S	EB8122C005S	5	EB8135C005S	EB8136C005S	



Tibial Insert



Quadrofit Standard Tibial Insert (Sterile)						
PE - UHMW Cross linked						
Ref. Number						Thickness
0	1	2	3	4	5	
EB8151U008S	EB8151U108S	EB8151U208S	EB8151U308S	EB8151U408S	EB8151U508S	8
EB8151U010S	EB8151U110S	EB8151U210S	EB8151U310S	EB8151U410S	EB8151U510S	10
EB8151U012S	EB8151U112S	EB8151U212S	EB8151U312S	EB8151U412S	EB8151U512S	12
EB8151U015S	EB8151U115S	EB8151U215S	EB8151U315S	EB8151U415S	EB8151U515S	15
EB8151U018S	EB8151U118S	EB8151U218S	EB8151U318S	EB8151U418S	EB8151U518S	18

Quadrofit Posterior Stabilized Tibial Insert (Sterile)

PE - UHMW Cross linked

Ref. Number						Thickness
0	1	2	3	4	5	
EB8155U008S	EB8155U108S	EB8155U208S	EB8155U308S	EB8155U408S	EB8155U508S	8
EB8155U010S	EB8155U110S	EB8155U210S	EB8155U310S	EB8155U410S	EB8155U510S	10
EB8155U012S	EB8155U112S	EB8155U212S	EB8155U312S	EB8155U412S	EB8155U512S	12
EB8155U015S	EB8155U115S	EB8155U215S	EB8155U315S	EB8155U415S	EB8155U515S	15
EB8155U018S	EB8155U118S	EB8155U218S	EB8155U318S	EB8155U418S	EB8155U518S	18



Quadrofit MB Standard Tibial Insert (Sterile)

PE - UHMW Cross linked

Ref. Number			Thickness
0	2	4	
EB8251U008S	EB8251U208S	EB8251U408S	8
EB8251U010S	EB8251U210S	EB8251U410S	10
EB8251U012S	EB8251U212S	EB8251U412S	12
EB8251U015S	EB8251U215S	EB8251U415S	15
EB8251U018S	EB8251U218S	EB8251U418S	18



Quadrofit MB Posterior Stabilized Tibial Insert (Sterile)

PE- UHMW Cross linked

Ref. Number			Thickness
0	2	4	
EB8255U008S	EB8255U208S	EB8255U408S	8
EB8255U010S	EB8255U210S	EB8255U410S	10
EB8255U012S	EB8255U212S	EB8255U412S	12
EB8255U015S	EB8255U215S	EB8255U415S	15
EB8255U018S	EB8255U218S	EB8255U418S	18



Tibial Component

Quadrofit Standard Tibial Component (Sterile)

Cemented CoCr	Size	Cementless CoCr Ti+Ha Coated
EB5361C000S	0	EB5362C000S
EB5361C001S	1	EB5362C001S
EB5361C002S	2	EB5362C002S
EB5361C003S	3	EB5362C003S
EB5361C004S	4	EB5362C004S
EB5361C005S	5	EB5362C005S



Quadrofit MB Tibial Component (Sterile)


Cemented CoCr	Size	Cementless CoCr Ti+Ha Coated
EB8261C000S	0	EB8262C000S
EB8261C001S	1	EB8262C001S
EB8261C002S	2	EB8262C002S
EB8261C003S	3	EB8262C003S
EB8261C004S	4	EB8262C004S
EB8261C005S	5	EB8262C005S



Tibial Component Half Spacer

Quadrofit Tibial Component Half Spacer (Sterile)

5mm Ti	Size	10mm Ti
EB5365T000S	0	EB5367T000S
EB5365T001S	1	EB5367T001S
EB5365T002S	2	EB5367T002S
EB5365T003S	3	EB5367T003S
EB5365T004S	4	EB5367T004S
EB5365T005S	5	EB5367T005S



Patella Component

Quadrofit Insert Patella Component (Sterile) PE - UHMW Cross linked	
Ref. Number	Size
EB8171U000S	0
EB8171U001S	1
EB8171U002S	2



Quadrofit ReSurfacing Patella Component (Sterile) PE - UHMW Cross linked	
Ref. Number	Size
EB8173U000S	0
EB8173U001S	1
EB8173U002S	2
EB8173U003S	3





Tibial Stem



Quadrofit Tibial Stem (Sterile)		
Ref. Number		Size
CoCr	Ti	
EB8291C025S	EB5391T025S	25
EB8291C050S	EB5391T050S	50
EB8291C080S	EB5391T080S	80
EB8291C100S	EB5391T100S	100

4+1 Quadrofit Total Knee System Instrument Set (EK1103N001N)



Tray 1



Tray 1

No	Ref.No	Description			
1	EE0851S205N	TIBIAL COMPONENT IMPACTOR			
2	EE0851S325N	TIBIAL STEM TEST FOR TOTAL KNEE SYSTEM / 25			
3	EE0852S033N	PUNCH HANDLE			
4	EE0851S350N	TIBIAL STEM TEST FOR TOTAL KNEE SYSTEM / 50			
5	EE0851N805N	POSTERIOR STABILIZATION ATTACHMENT FOR TIBIAL INSERT TEST			
6	EE0851S010N	TIBIAL COMPONENT TEST HANDLE FOR TOTAL KNEE SYSTEM			
7	EE0703S501N	GUIDE FOR D10 DRILL(BARREL NEW DESIGN)			
8	EE0852S045N	PUNCHER / 4 - 5			
9	EE0852S043N	PUNCHER / 2 - 3			
10	EE0852S041N	PUNCHER / 0 - 1			
11	EE0851S000N	TIBIAL COMPONENT TEST FOR TOTAL KNEE SYSTEM / 0			
12	EE0851S001N	TIBIAL COMPONENT TEST FOR TOTAL KNEE SYSTEM / 1			
13	EE0851S002N	TIBIAL COMPONENT TEST FOR TOTAL KNEE SYSTEM / 2			
14	EE0851S003N	TIBIAL COMPONENT TEST FOR TOTAL KNEE SYSTEM / 3			
15	EE0851S004N	TIBIAL COMPONENT TEST FOR TOTAL KNEE SYSTEM / 4			
16	EE0851S005N	TIBIAL COMPONENT TEST FOR TOTAL KNEE SYSTEM / 5			
17	EE0855N008N	TIBIAL INSERT TEST / 0 - 8	20	EE0855N308N	TIBIAL INSERT TEST / 3 - 8
	EE0855N010N	TIBIAL INSERT TEST / 0 - 10		EE0855N310N	TIBIAL INSERT TEST / 3 - 10
	EE0855N012N	TIBIAL INSERT TEST / 0 - 12		EE0855N312N	TIBIAL INSERT TEST / 3 - 12
	EE0855N015N	TIBIAL INSERT TEST / 0 - 15		EE0855N315N	TIBIAL INSERT TEST / 3 - 15
	EE0855N018N	TIBIAL INSERT TEST / 0 - 18		EE0855N318N	TIBIAL INSERT TEST / 3 - 18
18	EE0855N108N	TIBIAL INSERT TEST / 1 - 8	21	EE0855N408N	TIBIAL INSERT TEST / 4 - 8
	EE0855N110N	TIBIAL INSERT TEST / 1 - 10		EE0855N410N	TIBIAL INSERT TEST / 4 - 10
	EE0855N112N	TIBIAL INSERT TEST / 1 - 12		EE0855N412N	TIBIAL INSERT TEST / 4 - 12
	EE0855N115N	TIBIAL INSERT TEST / 1 - 15		EE0855N415N	TIBIAL INSERT TEST / 4 - 15
	EE0855N118N	TIBIAL INSERT TEST / 1 - 18		EE0855N418N	TIBIAL INSERT TEST / 4 - 18
19	EE0855N208N	TIBIAL INSERT TEST / 2 - 8	22	EE0855N508N	TIBIAL INSERT TEST / 5 - 8
	EE0855N210N	TIBIAL INSERT TEST / 2 - 10		EE0855N510N	TIBIAL INSERT TEST / 5 - 10
	EE0855N212N	TIBIAL INSERT TEST / 2 - 12		EE0855N512N	TIBIAL INSERT TEST / 5 - 12
	EE0855N215N	TIBIAL INSERT TEST / 2 - 15		EE0855N515N	TIBIAL INSERT TEST / 5 - 15
	EE0855N218N	TIBIAL INSERT TEST / 2 - 18		EE0855N518N	TIBIAL INSERT TEST / 5 - 18



Tray 2



Tray 2

No	Ref.No	Description
1	EE0851S173N	BLADE RUNNER
2	EE0802S082N	TIBIA INTRAMEDULLARY GUIDE
3	EE0802S084N	TIBIA EXTRAMEDULLARY GUIDE
4	EE0802S092N	TIBIAL ANKLE CLAMP
5	EE0802S090N	TIBIAL CUTTING BLOCK 5° - LEFT
6	EE0802S088N	TIBIAL CUTTING BLOCK 0° - LEFT
7	EE0802S089N	TIBIAL CUTTING BLOCK 5° - RIGHT
8	EE0802S087N	TIBIAL CUTTING BLOCK 0° - RIGHT
9	EE0802S085N	TIBIAL AXIS CONTROL PIECE
10	EE0802S086N	TIBIA STYLUS
12	EE0802S183N	TIBIA INTRAMEDULLARY GUIDE SUPPORT PIECE-ADAPTOR
13	EE0802S091N	TIBIAL DISTAL BODY
[14]	EE0401N114N	SCREWDRIVER Stryker 3,5 - optional
14	EE0401N065N	SCREW DRIVER / 3,5



Tray 3



Tray 3

No	Ref.No	Description
1	EE0715N005N	UNIVERSAL FEMORAL COMPONENT IMPACTOR
2	EE0852S032N	SLIDE HAMMER FOR TOTAL KNEE SYSTEM
3	EE0852S105N	FEMORAL COMPONENT IMPACTOR
4	EE0110S051N	NAIL HOLDER PLIERS
5	EE0651S005N	PIN ADAPTER FOR QUADROFIT REVISION KNEE SYSTEM
6	EE0651S071N	DRILL FOR TOTAL KNEE SYSTEM / Ø7 - Stryker Type
7	EE0651S101N	DRILL FOR TOTAL KNEE SYSTEM / Ø10 - Stryker Type
8	EE0802S056N	INTRAMEDULAR ROD FOR FEMUR
9	EE0802S055N	INTRAMEDULAR ROD FOR TIBIA
10	EE0856N010N	POSTERIOR STABILISATION SHAFT FOR FEMORAL COMPONENT TEST
11	EE0657N140N	NAIL FOR TOTAL KNEE SYSTEM / Ø 2,5 X 40 - 4pcs
12	EE0657N255N	NAIL FOR TOTAL KNEE SYSTEM / Ø3,0 X 55 - 6pcs
13	EE0658N001N	STEPPED NAIL / Ø 3,1 X 40 - 2pcs (New)
14	EE0658S001N	STEPPED THREADED NAIL / Ø 3,1 X 40 - 2pcs (New)
15	EE0652N101N	THREADED NAIL FOR TOTAL KNEE SYSTEM / Ø 3,2 X 70 - 4pcs
16	EE0657N570N	NAIL FOR TOTAL KNEE SYSTEM / Ø 3,2 X 70 - 4pcs



Tray 4



Tray 4

No	Ref.No	Description
1	EE0856N100N	FEMORAL COMPONENT TEST / RIGHT - 0
2	EE0856N101N	FEMORAL COMPONENT TEST / RIGHT - 1
3	EE0856N102N	FEMORAL COMPONENT TEST / RIGHT - 2
4	EE0856N103N	FEMORAL COMPONENT TEST / RIGHT - 3
5	EE0856N104N	FEMORAL COMPONENT TEST / RIGHT - 4
6	EE0856N105N	FEMORAL COMPONENT TEST / RIGHT - 5
7	EE0856N200N	FEMORAL COMPONENT TEST / LEFT - 0
8	EE0856N201N	FEMORAL COMPONENT TEST / LEFT - 1
9	EE0856N202N	FEMORAL COMPONENT TEST / LEFT - 2
10	EE0856N203N	FEMORAL COMPONENT TEST / LEFT - 3
11	EE0856N204N	FEMORAL COMPONENT TEST / LEFT - 4
12	EE0856N205N	FEMORAL COMPONENT TEST / LEFT - 5
13	EE0802S079N	DISTAL FEMUR GUIDE FOR COMPACT KNEE SYSTEM / 9°
14	EE0802S076N	DISTAL FEMUR GUIDE FOR COMPACT KNEE SYSTEM / 6°
15	EE0802S075N	DISTAL FEMUR GUIDE FOR COMPACT KNEE SYSTEM / 5°
16	EE0802S073N	DISTAL FEMUR GUIDE FOR COMPACT KNEE SYSTEM / 3°
17	EE0802S070N	DISTAL FEMUR GUIDE FOR COMPACT KNEE SYSTEM / 0°
18	EE0852S013N	CUTTING BLOCK STABILIZATION ATTACHMENT FOR COMPACT KNEE SYSTEM - 2pcs
19	EE0142S001N	CUTTING BLOCK HOLDER FOR COMPACT KNEE SYSTEM
20	EE0858S005N	FEMORAL CUTTING BLOCK FIXING SCREW FOR COMPACT KNEE SYSTEM - 2pcs
21	EE0860S001N	FEMORAL DISTAL CUTTING GUIDE ATTACHMENT FOR COMPACT KNEE SYSTEM
22	EE1202N000N	FEMORAL CUTTING BLOCK T FOR COMPACT KNEE SYSTEM / 0
23	EE1202N001N	FEMORAL CUTTING BLOCK T FOR COMPACT KNEE SYSTEM / 1
24	EE1202N002N	FEMORAL CUTTING BLOCK T FOR COMPACT KNEE SYSTEM / 2
25	EE1202N003N	FEMORAL CUTTING BLOCK T FOR COMPACT KNEE SYSTEM / 3
26	EE1202N004N	FEMORAL CUTTING BLOCK T FOR COMPACT KNEE SYSTEM / 4
27	EE1202N005N	FEMORAL CUTTING BLOCK T FOR COMPACT KNEE SYSTEM / 5
28	EE0859S001N	FEMORAL DISTAL CUTTING BLOCK FOR COMPACT KNEE SYSTEM
29	EE0852S016N	ANTERIOR SIZE GAUGE FOR COMPACT KNEE SYSTEM



Tray 5



Tray 5

No	Ref.No	Description
1	EE0891S000N	BOX GUIDE FOR FEMORAL COMPONENT - 0
2	EE0891S001N	BOX GUIDE FOR FEMORAL COMPONENT - 1
3	EE0891S002N	BOX GUIDE FOR FEMORAL COMPONENT - 2
4	EE0891S003N	BOX GUIDE FOR FEMORAL COMPONENT - 3
5	EE0891S004N	BOX GUIDE FOR FEMORAL COMPONENT - 4
6	EE0891S005N	BOX GUIDE FOR FEMORAL COMPONENT - 5
7	EE0802S200N	KNEE GAP GAUGE / 0-08
8	EE0802S201N	KNEE GAP GAUGE / 1-10
9	EE0802S202N	KNEE GAP GAUGE / 2-12
10	EE0802S203N	KNEE GAP GAUGE / 3-15
11	EE0802S204N	KNEE GAP GAUGE / 3-18
12	EE1041S002N	FEMORAL LIGAMENT REAMER - Stryker Type
13	EE0852S004N	NOTCHING INSTRUMENT
14	EE0861S000N	MB TIBIAL COMPONENT TEST / 0
	EE0861S001N	MB TIBIAL COMPONENT TEST / 1
	EE0861S002N	MB TIBIAL COMPONENT TEST / 2
	EE0861S003N	MB TIBIAL COMPONENT TEST / 3
	EE0861S004N	MB TIBIAL COMPONENT TEST / 4
	EE0861S005N	MB TIBIAL COMPONENT TEST / 5
15	EE0862S033N	MB PUNCH HANDLE ATTACHMENT
16	EE0861N008N	MB TIBIAL INSERT TEST / 0 - 8
	EE0861N010N	MB TIBIAL INSERT TEST / 0 - 10
	EE0861N012N	MB TIBIAL INSERT TEST / 0 - 12
	EE0861N015N	MB TIBIAL INSERT TEST / 0 - 15
	EE0861N018N	MB TIBIAL INSERT TEST / 0 - 18
17	EE0861N208N	MB TIBIAL INSERT TEST / 2 - 8
	EE0861N210N	MB TIBIAL INSERT TEST / 2 - 10
	EE0861N212N	MB TIBIAL INSERT TEST / 2 - 12
	EE0861N215N	MB TIBIAL INSERT TEST / 2 - 15
	EE0861N218N	MB TIBIAL INSERT TEST / 2 - 18
18	EE0861N408N	MB TIBIAL INSERT TEST / 4 - 8
	EE0861N410N	MB TIBIAL INSERT TEST / 4 - 10
	EE0861N412N	MB TIBIAL INSERT TEST / 4 - 12
	EE0861N415N	MB TIBIAL INSERT TEST / 4 - 15
	EE0861N418N	MB TIBIAL INSERT TEST / 4 - 18

4+1 Quadrofit Total Knee Instrument Set_Patella (EK1103N001P)



PATELLA

No	Ref.No	Description
1	EE0853N540N	BARREL FOR PATELLA CLAMP / 0
2	EE0853N541N	BARREL FOR PATELLA CLAMP / 1
3	EE0852S105N	BARREL FOR PATELLA CLAMP / 2
4	EE0853N550N	PATELLA REAMER / 0
5	EE0853N551N	PATELLA REAMER / 1
6	EE0853N552N	PATELLA REAMER / 2
7	EE0853N507N	PATELLA CLAMP ATTACHMENT WITH POINTED JAW
8	EE0853N511N	PATELLA CUTTING GUIDE / RIGHT
9	EE0853N512N	PATELLA CUTTING GUIDE / LEFT
10	EE0853N525N	PATELLA PRESSURIZING JAW
11	EE0853N528N	DRILLING TEMPLATE
12	EE0853S005N	CALLIPER GAUGE
13	EE0651S050N	DRILL FOR TOTAL KNEE SYSTEM / Ø5
14	EE0853N561N	PATELLA STYLUS
15	EE0853N505N	PATELLA CLAMP
16	EE0712S001N	IMPACTOR PIM FOR TOTAL KNEE SYSTEM
17	EE0857N100N	PATELLA COMPONENT TEST / 0
18	EE0857N101N	PATELLA COMPONENT TEST / 1
19	EE0857N102N	PATELLA COMPONENT TEST / 2
20	EE0857N200N	RE-SURFACING PATELLA COMPONENT TEST / 0
21	EE0857N201N	RE-SURFACING PATELLA COMPONENT TEST / 1
22	EE0857N202N	RE-SURFACING PATELLA COMPONENT TEST / 2
23	EE0857N203N	RE-SURFACING PATELLA COMPONENT TEST / 3



HOW TO CONTACT US

Euromed Implants GmbH
Bei St. Wilhadi 1
21682 Stade, Germany

Phone: +49 (0) 4141 7897 290
Fax: +49 (0) 4141 7782 700

office@euromed-implants.de
www.euromed-implants.de